

PLAN EVALUATION IN IMRT

Prof G. KILARA
Head Dept Of Radiation Oncology
Curie centre of oncology
Bangalore

“A dose plan and treatment delivery that is optimized using inverse or forward planning techniques for modulated beam delivery, using EITHER a binary collimator OR with a conventional MLC system, using either “Sliding Window (DMLC) OR “Step and Shot” (SMLC) modes”

National Cancer Institute
Guidelines for IMRT in
Clinical Trials

AIM OF RADIATION THERAPY

- TUMORICIDAL DOSE TO TUMOUR
- LEAST DOSE TO SURROUNDING
NORMAL TISSUE

CONVENTIONAL RADIOTHERAPY

- RECTANGULAR / SQUARE FIELDS
- LIMITED NO. OF FIELDS
- UNIFORM DOSE INTENSITIES IN THE
FIELD

CONVENTIONAL RADIOTHERAPY

NORMAL TISSUE SPARING :

- WEDGE FILTERS
- BEAM MODIFICATION DEVICES
- SHIELDING BLOCKS
- X RAY SIMULATION
- IMMOBILISATION etc...

CRITICAL PARAMETERS

- PRECISION IN DOSE DELIVERY
- REPRODUCIBILITY
- TUMOUR AND NORMAL TISSUE
DELINEATION

ADVANCES


- 3D RECONSTRUCTIONS
(CT/MRI/PET)
- 3D COMPUTERISED TREATMENT
PLANNING SYSTEM
- MULTILEAF / MICROMULTILEAF
COLLIMATION

ABILITY NOW TO

- VISUALISE STRUCTURES IN 3D
- SHAPE IRREGULAR RADIATION FIELDS
IN 3D
- MATCH THE VISUALISED TUMOUR
SHAPE AND RADIATION DOSE SHAPE
TO EACH OTHER

LATEST TOOL TO ENHANCE THE MATCH
BETWEEN TUMOUR AND DOSE SHAPE AND
SPARE NORMAL TISSUE:

INTENSITY MODULATED
RADIATION
THERAPY (IMRT)



INTENSITY MODULATED RADIATION THERAPY (IMRT) REQUIREMENTS

- **LINEAR ACCELERATOR 6 MV**
- **MULTILEAF/MICROMULTILEAF
COLLIMATOR**
- **3D TPS WITH CT/ MRI / PET FUSION**
- **INVERSE – PLANNING SOFTWARE**
- **IMMOBILISATION**
- **SPECIALIZED IMRT QA EQUIPMENT**

INVERSE PLANNING

COMPUTER ALGORITHM ALLOWS TO SPECIFY
DOSE CONSTRAINTS TO NORMAL TISSUE (eg
BLADDER, RECTUM, CORD etc) WHILE GIVING
TUMOUR REQUIRED DOSE


[Cf Forward planning]

IMRT INDICATIONS


- AS SOLE MODALITY
- AS BOOST TREATMENT

(PROSTATE, HEAD AND
NECK, BRAIN, BREAST, CERVIX,
PARASPINAL TUMOURS , LUNG,
BLADDER etc..)

PROSTATE CANCER




PROSTATE CANCER



Skin
bones


3D RECONSTRUCTION

PROSTATE CANCER



3D RECONSTRUCTION

PROSTATE CANCER



Skin
bones
intestine
rectum
Bladder
prostate

TREATMENT BEAMS


Arc No. 7 Step No. 1

PROSTATE CANCER

Gantry

Skin
bones
intestine
rectum
Bladder
prostate

MICRO MULTI LEAF COLLIMATION




0010_2005, Treatment Name: IMAT_7
Zoom: 0.5

Isocenter


Xf = -48.7

Yf = -4.9


Zf = -21.3



Dose Volume Histogram - IMAT_3



Dose Volume Histogram



Hide Toolbar

Hide Legend

Object Name

rectum

Dose relative to:

- ☐ Treatment Maximum
- ☒ Absolute
- ☐ Prescribed

Name	Show
conformal_3d	Show
IMAT11	Hide
stepandshoot	Hide
IMAT_c	Hide
stepandshoot_c	Hide
IMRT001a	Hide
IMAt_1	Hide

MaxDose [cGy]

1939

Print

Print

- ☒ To Printer
- ☐ To File

IMRT PLAN EVALUATION PARAMETERS

A . CLINICAL PARAMETERS

- Patient selection
- Site, Stage
- Psychology
- 3DCRT Vs IMRT
- Target delineation
- Immobilization
- Dose prescription guidelines

IMRT PLAN EVALUATION PARAMETERS

B . PHYSICAL PARAMETERS

- TPS (Algorithms etc.)
- Dose constraints
- DVH
- GTV $\longrightarrow$ PTV (setup errors)
- OARs
- QA

IMRT PLAN EVALUATION PARAMETERS

C. EQUIPMENT PARAMETERS

- Linear Accelerator
- Type of MLC / μ MLC
- Mould room
- Simulation
- Port films
- MU (DR)
- TPS

CLINICAL CONSIDERATIONS

- IMMOBILIZABLE SITE
- IRREGULARLY SHAPED CONCAVE EDGES
- ADJACENT TO CRITICAL ORGANS
- COOPERATIVE PATIENTS

IDEAL FOR
IMRT

Morbidity obese, poor immobilization
Claustrophobic patient, moribund

POOR
CANDIDATES

CLINICAL CONSIDERATIONS

“HOT” and “COLD” SPOTS IN 3D.

- MAGNITUDE
- VOLUME
- LOCATION

“COLD” - NOT WITHIN GTV - AS FAR AWAY.

“HOT” - DESIRED IN-HOMOGENITY

CLINICAL CONSIDERATIONS

FOLLOW GUIDELINES : Egs :

HEAD & NECK : RTOG H-022 TRIAL
95% PTV BY PRESCRIPTION
ISODOSE

GYNAEC : RTOG 0418

PHYSICS PARAMETERS

IMRT vs 3DCRT

- LARGE VOLUME OF NORMAL TISSUE RECEIVE LOW DOSES
- ↑ MONITORING UNITS NEEDED FOR DELIVERY
- VERIFICATION PORT FILMS REQUIRED
- ↑ PATIENT EXPOSURE FROM HEAD LEAKAGE, NEUTRON PRODUCTION AND SCATTER

? ↑ RISK OF 2nd MALIGNANCY

?? ↑ RT PNEUMONITIS IN THORACIC IMRT

PHYSICS PARAMETERS

- BETTER FOR CONCAVE SURFACE S
- PLANNING IS LESS DEPENDENT ON BEAM ENERGY
- CONFORMITY AT THE EXPENSE OF LOW DOSE TO LARGE VOLUME.

ie CHECK 30% ISODOSE

PHYSICS PARAMETERS

- PLANNING IS TIME CONSUMING
ie NOT FOR RAPID PLANS
Eg . EMERGENCY Rx etc.
- ALGORITHM FOR ITERATION IMPORTANT
HETEROGENITY CORRECTION TO BE USED

PHYSICS PARAMETERS

- AVOID MATCHING TWO IMRT FIELDS OR
IMRT +PHOTONS/ELECTRON
- ORGAN MOTION AND SET UP ERRORS CAUSE
DOSIMETRIC PROBLEMS
(PATIENTS MOVES / ORAGAN MOVES)

IMRT PROSTATE DOSE VOLUME CONSTRAINTS – OAR

AUTHOR	RECTUM	BLADDER	OTHERS
1.Zelefsky et al	$\leq 30\%$ to get ≥ 75.6 $\leq 53\%$ to get ≥ 47	$\leq 53\%$ to get ≥ 47	NS
2.Ezzell et al	$\leq 40\%$ to get ≥ 65 $\leq 30\%$ to get ≥ 70 $\leq 10\%$ to get ≥ 75	$\leq 30\%$ to get ≥ 75 $D_{max} \leq 81$	FH $D_{max} \leq 50$
3. Sethil et al	$\leq 30\%$ to get ≥ 65	$\leq 30\%$ to get ≥ 65	PP $\leq 25\%$ To get $\geq 40\%$

IMRT PROSTATE TARGET DEFINITIONS

	CTV	PTV	Prescription (TD/ Fr) Gy
1. Zelefsky 81 Gy plan 86.4 Gy plan	P+ SV P + SV	CTV +1 cm UE CTV + 1 cm	PTV 81/1.8 ≥90% to get ≥70 PTV 86.4/1.8 ≥85% to get ≥86.4
2. Ezzell	P + SV	CTV + 1 cm UE	75.6/1.8 to get ≥95% CTV
3. Sethil et al		PTV1=(P +SV)+1 cm UE PTV2 =(p) + 1 cm UE	PTV1 55.8/1.8 PTV2 18/1.8 ; 25.2/1.8

EQUIPMENT PARAMETERS

- BEAM DIRECTION IMPORTANT
- AVOID TREATMENT THROUGH IMMOBILIZATION DEVICES WHICH ATTENUATE BEAM
- AVOID LONG TREATMENT TIME
 - INCOVENIENT
 - ? RADIOBIOLOGICAL
- USE APPROPRIATE TPS

EQUIPMENT PARAMETERS

- STRINGENT QA IS CENTRAL TO SUCCESS OF IMRT
 - PROPER IMMOBILIZATION IS CRUCIAL
 - LASERS IN CT SIMULATOR IS MUST
 - FLAT CT COUCH FOR SIMULATION
- SET UP ERRORS – SYSTEMATIC & RANDOM
- QA OF LINAC AND TPS IS ESSENTIAL

FINALLY

IMRT IS AN EXCELLENT TECHNIQUE IF USED
APPROPRIATELY BY A COMPETENT AND EXPERIENCED
TEAM USING THE NECESSARY HIGH QUALITY MOULD
ROOM AND SIMULATION PARAMETERS AND IMPLIMENTED
UNDER A STRINGENT QA PROGRAMME

REMEMBER : G – 1- G.0 !

THANK YOU....